
WIDE RANGE OF MOTION ADJUSTMENT

Brace is shipped with 10° hinge extension stops installed with additional in and extension stops in place. Stops are labeled in degrees:

Extension stops: 0°, 10°, 20°, 30° & 40°
Flexion stops: 45°, 60°, 75° & 90°

Brace is shipped with no flexion stops installed

To change hinge stops, first remove the condyle pads held in place by hook/loop fasteners. Simply pull the pads off the hinge plates. To replace the pads be that they are positioned properly and firmly push them onto their hook/loop fasteners.

Use a screwdriver to remove screws. Remove existing stop, insert new replace screw. Be sure that screws are seated and snug. Both medial (left) hinge and lateral (outside) hinge stops must be identical in degrees.

OSTEOARTHRITIS (OA) HINGE INSTRUCTIONS

With the brace properly applied and tightened, OA pain symptoms can be relieved by applying a lateral force to the outside of the knee by turning the dial (in the hinge) clockwise. With the knee bent in either a seated or lying position, you can easily turn the circular dial to move the condyle pad (toward knee). The window on the top edge of the dial will display gings ranging from 0 to 8 in 0.5mm increments; the higher the number the greater the lateral force being applied. To reduce the lateral force, turn the dial counterclockwise to a lower setting.

Make sure the large gel condyle pad is attached to the lateral (outside) condyle basket.

CARE AND CLEANING

Brace must be removed and cleaned regularly depending upon your activity and perspiration. Hand wash straps and frame pads in cold water with detergent, rinse, AIR DRY. Occasionally apply a dry lubricant to hinge daily after use in water.

Frame pads are held in place by hook/loop fasteners. Carefully pull the off the brace frame. To replace the pads be sure that they are positioned properly and firmly push them onto their hook/loop fasteners.

Condyle hinge condyle pads are held in place by hook/loop fasteners. Simply pull pads off the hinge plates. To replace the pads be sure that they are seated properly and firmly push them onto their hook/loop fasteners.

To remove the hook/loop straps first remove the strap pads on the inside of each. Then loosen the hook/loop tabs at each strap end and feed the straps through the anchor points on the side of the brace until the strap is free. Be sure the strap pads and strap pads according to their strap number locations, illustration.

TIGHTENING

To achieve the best possible fit of the brace straps, be sure that the strap pads inside of each strap are centered on each strap and provide a comfort contact point between the skin and the strap.

Length Adjustment: If a strap is too long or too short (preventing the strap from properly adhering to the strap). First, remove the strap pad and unfasten the strap tab, then unfasten the "tail end" of the strap slide the strap to either lengthen or shorten. Refasten the strap ends and use the strap pads.

If no diagnosed or apparent instability associated with the unicompartmental knee OA, simply the strap application for FUSION OA braces by using all straps behind the knee first (1, 2, 3, and 5) and then apply straps 4 in front of the knee.

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